

BI@BM1000

PORTABLE DUAL-HEAD VEIN FINDER



www.portable-ultrasound-scanner.com

Quality, Professional, Innovational

BM1000

PORTABLE DUAL-HEAD VEIN FINDER

Product Appearance



Product Features

- Using infrared light source, no radiation, no laser
- Flexible screen with high resolution, can be placed to the most appropriate position
- Better penetration, working for hand, arm, elbow, foot and head
- Focus adjustable, easy to get clear image
- Bone bracket capable of bending of arbitrary shape, angle
- Internal lithium battery, can supply power for more than 1 hour
- Exquisite suitcase, protective and portable, as standard configuration


Product Applications

- Hospital, pediatric clinic, injection room
- Outpatient clinic, medical laboratory
- Blood bank, military medics unit, the CDC
- Rehabilitation center, nursing homes and others


Product Specifications

- Power supply: AC100-240V 50/60Hz 15VA
- Input power: DC 12.6V 1A
- Internal power supply: Lithium battery DC12V 3000MAH
- Voltage: DC11.1V-12.6V
- Inspecting distance: 10cm-30cm(adjustable)

Actual Images




leg



dorsal feet



dorsal hand



dorsal hand

- Screen size: 5 inches(resolution:800*480)
- Net weight: 3.5Kg
- Gross weight: 4Kg
- Dimension (package): 43 X 27 X 25cm

Wuxi Biomedical Technology Co., Ltd. www.portable-ultrasound-scanner.com

Add:Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel:86-510-85572789

Fax:86-510-85572789

Email:sales@cnbiomedical.com

Mobilephone:86-13382885910

BIC is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.