

Projection Type Vein Finder

Product Feature

- Near infrared light vein finder, no laser, no radiation.
- High contrast image, no noise or interference.
- Vein image projected to human body directly.
- No skin touch: real-time image on skin surface.
- Easy to operate: one button to power on/off.
- Easy to carry: smart and flexible.
- Built-in rechargeable lithium battery: Continuously work for 4 hours without power supply.
- Adapter to AC/DC power.
- Projection color, size and brightness are adjustable.

Specification

- Light type : Near infrared light
- Visible vein size : $\geq 1\text{mm}$
- Detecting depth : 10mm
- Accuracy : $\leq 0.3\text{mm}$
- Best imaging distance : 20-26cm
- Charging time : 2 hours
- Appearance : blue, gray, pink
- Handheld/Table/Trolley type

www.portable-ultrasound-scanner.com

Quality, Professional, Innovational

Product Picture





Product Applications

- Hospital, pediatric clinic, injection room
- Outpatient clinic, medical laboratory
- Blood bank, military medics unit, the CDC
- Rehabilitation center, nursing homes and others

Application Scope

- Children
- Infants
- The elderly
- Obesity
- Dark skin
- Patients with anemia, hypotension, or excessive loss of blood
- Plastic surgery
- Vascular surgery



Wuxi Biomedical Technology Co., Ltd.

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Email: sales@cnbiomedical.com

Fax: 86-510-85572789

Mobilephone: 86-13382885910

www.portable-ultrasound-scanner.com



is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.