

BIO+ HSL-NEW
Digital Hand-held Slit Lamp



www.portable-ultrasound-scanner

Quality, Professional, Innovational

Feature And Advantages

- 3.97 inch TFT-LCD color touch screen
- Real-time image display, 10 times magnification
- Brand new optical system, larger field of view
- Video function offering live data
- Easy to operate and high definition imaging
- Multi data transfer and storage system
- Light weight 850g, movable and portable
- Battery works for 6 hours continuously



3 storage modes




- Built-in 16G storage(can storage 80,000 images)
- 2USB connect the computer and export.
- Real-time Wifi, computer, mobile terminal and cloud transmission, as well as data analysis and software upgrading services.

Product Parameter

Model	HSL-NEW
Screen	3.97 inch TFT-LCD color touch screen
Illumination angle	≥30°
Magnification	10X
Working distance	80mm
Slit width	0mm~12mm continuously adjustable
Object eyepiece	0.42X
Aperture diameter	0.2mm,1mm, 5mm, 12mm
Filter	Heat-absorption, Red-free, Cobalt blue
Light source	White LED
Image resolution	16MP
Transmission	Mini USB,WIFI,BT
Storage	16GB MicroSD card
Power	3.6V/3400mAh rechargeable battery
Forehead rest bar	≥10mm
Numerical aperture	0.05
Interface	Micro USB
Focus Mode	Auto
N. Weight	850g(Typical)

Real Images



Wuxi Biomedical Technology Co., Ltd.

www.portable-ultrasound-scanner.com

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: sales@cncbiomedical.com

Mobile phone: 86-13382885910

BIO is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.